

T.A. 01

Multi-parameter titration unit

This titration unit is designed to make easier the determinations by adding a titration solution with a changing colored indicator.

The unit is composed of a 10 ml automatic zero setting burette in 1/20 ml, managed by an electrovalve, a control station box including a magnetic agitator with variable speed, and a retro lighted tray which implies a better visualization of the end of the titration.

Advantages

- Rapid starting
- Easy to use
- Excellent visualization of the end of the titration thanks to the retro lighted tray
- Precision and reliability
- Very competitive price

Specific features

- Dimensions
 - Height : 670 mm
 - Width : 185 mm
 - Depth : 230 mm
- Weight : 1,9 kg
- Power supply: 220 V - 50 Hz
- 10 ml burette in 1/20 ml
- Flow control by clip-on electrovalve
- Automatic zero setting
- Magnetic agitator with variable speed
- Retro lighted tray



Ref. 114000

LABORATOIRES DUJARDIN-SALLERON

872, route de la Gare
 37210 NOIZAY - FRANCE
 Tél : +33 (0)2 47 25 58 25
 Fax : +33 (0)2 47 25 58 30
 info@dujardin-salleron.com
 www.dujardin-salleron.com